

NAME

DBI::Gofer::Serializer::Storable - Gofer serialization using Storable

SYNOPSIS

```
$serializer = DBI::Gofer::Serializer::Storable->new();

$string = $serializer->serialize( $data );
($string, $deserializer_class) = $serializer->serialize( $data );

$data = $serializer->deserialize( $string );
```

DESCRIPTION

Uses *Storable::nfreeze()* to serialize and *Storable::thaw()* to deserialize.

The *serialize()* method sets local `$Storable::forgive_me = 1`; so it doesn't croak if it encounters any data types that can't be serialized, such as code refs.

See also `DBI::Gofer::Serializer::Base`.