

NAME

EC_GFp_simple_method, EC_GFp_mont_method, EC_GFp_nist_method,
 EC_GFp_nistp224_method, EC_GFp_nistp256_method, EC_GFp_nistp521_method,
 EC_GF2m_simple_method, EC_METHOD_get_field_type - Functions for obtaining
 EC_METHOD objects.

SYNOPSIS

```
#include <openssl/ec.h>

const EC_METHOD *EC_GFp_simple_method(void);
const EC_METHOD *EC_GFp_mont_method(void);
const EC_METHOD *EC_GFp_nist_method(void);
const EC_METHOD *EC_GFp_nistp224_method(void);
const EC_METHOD *EC_GFp_nistp256_method(void);
const EC_METHOD *EC_GFp_nistp521_method(void);

const EC_METHOD *EC_GF2m_simple_method(void);

int EC_METHOD_get_field_type(const EC_METHOD *meth);
```

DESCRIPTION

The Elliptic Curve library provides a number of different implementations through a single common interface. When constructing a curve using `EC_GROUP_new` (see [EC_GROUP_new\(3\)](#)) an implementation method must be provided. The functions described here all return a const pointer to an `EC_METHOD` structure that can be passed to `EC_GROUP_NEW`. It is important that the correct implementation type for the form of curve selected is used.

For F_2^m curves there is only one implementation choice, i.e. `EC_GF2_simple_method`.

For F_p curves the lowest common denominator implementation is the `EC_GFp_simple_method` implementation. All other implementations are based on this one. `EC_GFp_mont_method` builds on `EC_GFp_simple_method` but adds the use of montgomery multiplication (see [BN_mod_mul_montgomery\(3\)](#)). `EC_GFp_nist_method` offers an implementation optimised for use with NIST recommended curves (NIST curves are available through `EC_GROUP_new_by_curve_name` as described in [EC_GROUP_new\(3\)](#)).

The functions `EC_GFp_nistp224_method`, `EC_GFp_nistp256_method` and `EC_GFp_nistp521_method` offer 64 bit optimised implementations for the NIST P224, P256 and P521 curves respectively. Note, however, that these implementations are not available on all platforms.

`EC_METHOD_get_field_type` identifies what type of field the `EC_METHOD` structure supports, which will be either F_2^m or F_p . If the field type is F_p then the value `NID_X9_62_prime_field` is returned. If the field type is F_2^m then the value `NID_X9_62_characteristic_two_field` is returned. These values are defined in the `obj_mac.h` header file.

RETURN VALUES

All `EC_GFp*` functions and `EC_GF2m_simple_method` always return a const pointer to an `EC_METHOD` structure.

`EC_METHOD_get_field_type` returns an integer that identifies the type of field the `EC_METHOD` structure supports.

SEE ALSO

[crypto\(3\)](#), [ec\(3\)](#), [EC_GROUP_new\(3\)](#), [EC_GROUP_copy\(3\)](#), [EC_POINT_new\(3\)](#),
[EC_POINT_add\(3\)](#), [EC_KEY_new\(3\)](#), [d2i_ECPKParameters\(3\)](#), [BN_mod_mul_montgomery\(3\)](#)